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Lighting and Imaging
SCHOTT AG
 Otto-Schott-Straße 2
 55127 Mainz
 Germany
 Phone: +49(0)6131/66-7844
 Fax: +49(0)6131/66-7850
www.schott.com/lightingimaging

SCHOTT
 glass made of ideas

Universal Light Guide

Short Manual

Intended use:

The light guide is used to transmit light during endoscopic diagnosis and therapy or for other medical illumination applications when the properties and operating instructions are complied with. The light conductors are designed for use with halogen, xenon or LED based sources of cold light which are utilised in medical applications that involve endoscopes, medical instruments or microscopes and are equipped with an optical injection element for light conductors with a bundle diameter of 4.8 mm to 5.0 mm or 3.5 mm. High intensity visible light is conducted to the instrument through flexible light guides. The use of the light conductor together with LASER light sources and/or HF devices is not permitted.

Contraindications related to the product are currently not known.

Treatment:

Before first use and every subsequent use, the light conductor must be cleaned, disinfected, and/or sterilised by means of validated processes. Preparation is not permissible with the adapter sleeves screwed on. They must be removed and prepared separately prior to the preparation of the light conductor. The general recommendations / instructions on how to prepare medical products must be observed.

Manual procedures:

Coarse fouling must be removed beforehand as necessary. The light conductor must be placed in a suitable cleaning and disinfectant solution that covers all of its surfaces. The instructions of the disinfectant producers relating to concentration and residence time must be observed. Any adhering dirt must be wiped off when disinfection is complete. Following disinfection and cleaning, the light conductors must be thoroughly rinsed with fully desalinated water in order to remove any residue. Afterwards, use a clean, lint free cloth to dry the light conductors. Metal brushes, steel wool or abrasive cloths or sponges are not suitable for cleaning.

Mechanical procedures:

Use only validated preparation procedures whose efficiency has been proven for mechanical cleaning and disinfecting. The used cleaning agents and equipment (disinfectors) must be suitable for the intended purpose.

Steam sterilisation:

Steam sterilisation must be carried out in devices that perform validated procedures. Suitable procedures are the fractionated vacuum procedure and the gravitation procedure. Suitable procedures are performed at 134° C / minimum residence time 10 min. or 121° C / minimum residence time 20 min.

Low-temperature plasma (LTP) sterilisation:

- LTP steriliser in accordance with DIN EN ISO 14937
- validated according to DIN EN ISO 14937 (valid IQ/OQ (commissioning) and product specific performance qualification (PQ))
- maximum sterilisation temperature 60 °C (140 °F; including permissible tolerance)
- Long cycle (Sterrad 100S)

Reusability:

If handled with the necessary care and remaining undamaged and unsoiled, the light conductors can be reused up to 100 times; using them beyond this limit or using damaged and/or soiled instruments is the responsibility of the user. No liability in case of non observance.

The fact that the light conductor is reusable up to 100 times is validated in the following conditions:

Mechanical cleaning and disinfection in MIELE Professional PG8535 disinfectant using NOEDISHER MediClean forte (0.5 vol. %) cleaner; thermal disinfection at 95° C with a residence time of 10 min.

Sterilisation using a (triple) fractionated vacuum procedure in a LAUTENSCHLÄGER PROTOCERT 839 steam steriliser (134° C, 10 min., 3 bar).

For more information on how to use and prepare the light guide, detailed instructions are available for download at

www.schott.com/endoscopy-light-guide

Warning:

- Immediately before and after every use, the light conductor and accessories used with it must be examined for fault free state, correct functioning, and completeness.
- Damaged light conductors must be withdrawn from use immediately.
- Because of the high amounts of light energy, the light input surfaces and illuminated surfaces can heat up. Danger of burns!
- Light sources may be used for the above purpose only if the maximum permitted temperature of use of 350° C on the input side is complied with.
- Never look into the open end of a light conductor which is connected to a light source! Risk of blinding!
- The light conductor is not intended for invasive use.
- The light conductor is not provided in a sterile state.

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