Photo Foto: SCHOTT

“Double chamber cartridges offer more safety and convenience for the patient.”

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Andreas Wesp, SCHOTT Product manager

# A Double Chamber with a Fourfold Advantage

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Modern drugs often consist of combined active ingredients. Now these can be stored separately in a double chamber cartridge and yet be administered very easily.

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ALEXANDER LOPEZ

**F**or the organizers of Good Days, July 10 could be their day. This non-profit association from the United States is committed to proclaiming a “Chronic Disease Awareness Day,” an anniversary for all those who are permanently dependent on medications due to chronic diseases. “For chronic patients, easy access to

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Photo Foto : SCHOTT

**Active despite a chronic illness – modern pharmaceutical agents make this possible. Many of these agents must be liquefied prior to injection. SCHOTT has developed a solution for easier handling: innovative double chamber cartridges (right).**

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medication means being able to lead an independent life,” the organization notes. The pharmaceutical industry has made great progress in developing increasingly effective medicines for the chronically ill. These usually have to be injected, some on a weekly basis, others several times a day. That places a burden on patients as well as healthcare systems. It therefore makes sense to have the patients perform the injections themselves. This in itself is not a problem; however, many new drugs can only be stored stably in freeze-dried form. Prior to injection, they must be liquefied by using water for injection (WFI) or another diluent. The patient needs two vials with different contents and a syringe. First, the diluent must be drawn into the syringe and injected into the second vial. After shaking it, the drug can be drawn into the syringe again and administered. The risk of error, such as a possible contamination of the substances, increases with the number of

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steps. In addition, there is the risk that the patient could miss the right mixing ratio or inject the wrong dosage.

SCHOTT has developed a solution for these forms of administration: a double chamber cartridge for pen systems. It contains the two drug components in two consecutive chambers that are separated by a plunger. When the pen lock is turned, the plunger will be pushed to the level of the bypass, and the liquid will flow into the anterior chamber through the nearly 0.1 millimeter narrow opening. It will then mix with the second component. The final drug can then be injected without any risk of contamination and in the exact dosage by simply triggering the pen.

SCHOTT cartridges are compatible with the most commonly used pen systems and can be adapted to meet the specific needs of pharma companies. The position of the bypass can be changed or the edges can be modified, for example. "This is important if a customer wants to fill the drug directly in powder form," says SCHOTT Product Manager Andrea Wesp. Alternatively, the active ingredient can also be filled in liquid form and then be lyophilized inside a freeze dryer. The system is convenient for the patient and of great benefit to pharmaceutical manufacturers. "Many new drugs are initially made available as lyophilisates because manufacturers are still looking for a way to store the active ingredients in a stable manner," says Wesp. "If they opt for a double chamber, they don't have to accept any compromises because the components can be stored and protected, reconstituted quite easily and then administered in the safest and most convenient way – four advantages at once."

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
Source Quelle : schott/C3

SCHOTT double chamber cartridges are used in pens and can combine either liquids with liquids or liquids with medications in powder form. The two chambers are separated by a plunger. When the patient or doctor turns the closure of the pen, the back plunger moves forward and presses the medication against the front plunger up into a passageway that is less than 0.1 millimeters in width. The medications stored in the back chamber can pass through this bypass and enter the front chamber where they can mix. To inject the medication, the patient simply places the pen against his skin and presses the release button.

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