**DURAN-  
Laborglas  
Ausgabe I/2000**

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**DURAN  
Laboratory glassware  
Edition I/2000**

**All items are subject to technical  
and standards related changes.**

In cases of changes in DIN specifications, due, for example, to assimilation to ISO recommendations, our dimensions will be correspondingly adjusted over a reasonable time period.

Delivery through authorized dealers.

We feel unable to grant any rebates on orders of a net value below DM 500,-.

Articles for which standard pack quantities are given, can be supplied only in the quantities mentioned,

these quantities being at the same time the minimum number of pieces to be purchased. Exceptions to the above are specially marked.

Please make your arrangements accordingly.

The terms of delivery and payment are always applicable in their latest edition.

Reproduction of this catalogue, in whole or in part, is permitted only with our written consent.

As from the date this catalogue takes effect, all issues published previously will be void.

**Trademarks**

SCHOTT, DURAN, DUROPLAN, FIOLAX, KPG, SUPREMAX, PRODURAN, SBW, STERIPLAN,

AR-GLASS (soda-lime), KECK and WERTEX are registered trademarks of the SCHOTT-Group.

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## DURAN by Schott

DURAN is a borosilicate glass type 3.3 as specified in international standard ISO 3585 and DURAN products meet the most important international standards, for example the German<sup>1</sup>, British, American and French standards. Maximum chemical resistance, minimum thermal expansion and the resulting high resistance to thermal shock are its principal properties. This optimum physical and chemical

performance makes DURAN the ideal material for use in the laboratory and for large-scale industrial plants in the manufacture of chemical plant. Moreover, it is widely used on an industrial scale in all other areas of application in which extreme heat resistance, resistance to thermal shock, mechanical strength and exceptional chemical resistance are required.

### Definition

Glass is an inorganic mixture fused at high temperature which solidifies on cooling but does not crystallize. Its basic components, network formers and modifiers are present in the form of oxides in ordinary glasses. Typical glass formers (network formers) are silica ( $\text{SiO}_2$ ), boric acid ( $\text{B}_2\text{O}_3$ ), phosphoric acid ( $\text{P}_2\text{O}_5$ ) and, under certain circumstances, also aluminium oxide ( $\text{Al}_2\text{O}_3$ ).

These substances are capable of absorbing (dissolving) certain amounts of metal oxides without losing their glassy character. This means that the incorporated oxides do not participate as glass formers but modify certain physical properties of the glass structure as „network modifiers“.

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<sup>1</sup> DIN and EN standards can be purchased from: Beuth Verlag GmbH, Postfach 11 45, D-10772 Berlin  
ISO denomination according to ISO 4791/1

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## Glass structure

A large number of chemical substances solidify from the liquid state of the melt in the form of a glass. The glass formation depends on the cooling rate and requires mixed types of bonding (covalent bond and ionic bond) between the atoms or atom groups. As a result, glass-forming products have a strong tendency while still in the melt for three-dimensional crosslinking by polymerisation in a largely random manner. Crystals are formed when the individual atoms form a regular three-dimensional arrangement in what is known as a

“crystal lattice”, as soon as the particular substance changes from the liquid to the solid state. Glass, however, on cooling from the liquid state, forms a largely spatially random “network”. The main components which participate in the glass formation are therefore called “network formers”. Ions can be incorporated in this network of glass-forming molecules, as a result of which they tear up the network in certain places and modify the network structure and thus glass properties. That is why they are called “network modifiers”.

## Chemical composition of DURAN®

DURAN glass which is used in the laboratory and in industrial scale chemical plant, due to its excellent chemical and physical properties, has the following typical composition:

81 % by weight  $\text{SiO}_2$   
13 % by weight  $\text{B}_2\text{O}_3$   
4 % by weight  $\text{Na}_2\text{O}/\text{K}_2\text{O}$   
2 % by weight  $\text{Al}_2\text{O}_3$

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## Chemical properties

The chemical resistance of glass is more extensive than that of all other known materials. Borosilicate glass is highly resistant to water, acids, salt solutions, organic substances and also halogens such as chlorine and bromine. It also has relatively good resistance to alkaline solutions. Only hydrofluoric acid, concentrated phosphoric acid and strong alkaline solutions cause appreciable surface removal of the glass at higher temperatures.

### Hydrolytic resistance according to ISO 719 (98 °C)

Hydrolytic resistance by the glass grain method class ISO 719-HGB 1 (corresponds to the previous DIN 12111, hydrolytic class 1).

### Hydrolytic resistance according to ISO 720 (121 °C)

Hydrolytic resistance by the glass grain method class ISO 720-HGA 1.

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Fig. 2:  
Hydrolytic attack on DURAN as a  
function of time (100 °C) approx. to  
ISO 719.

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**Acid resistance according to DIN 12116**

Acid class 1

**Acid resistance according to DIN-ISO 1776**

Removal of  $\text{Na}_2\text{O} < 100 \mu\text{g}/\text{dm}^2$   
 The attack curves in Figure 3 show a maximum in acid ranges between 4-7 N (HCl, for example in the azeotrope containing 20.2% by weight). At higher concentrations, the reaction rate decreases significantly, so that the layer thicknesses which are

attacked are only in the range of a few thousandths of a millimetre after years. Thus, the mechanisms of acid attack are not relevant given the wall thicknesses used in practice. Therefore, it is rightly called "acid-resistant borosilicate glass".

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Fig. 3: Acid attack on borosilicate glass as a function of acid concentration.

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**Alkali resistance according to DIN 52322**

Alkali resistance class DIN 52322-A2 (corresponds to alkali resistance class ISO 695-A2). The surface removal curve in Figure 4 shows that visible attack on the glass surface does not take place until above about 60 °C. At lower temperatures (incorporation of coolers in the NaOH cycle), the reaction rates are so low that hardly any reduction of the wall

thickness takes place over a period of years. Long-term tests have shown that the use of NaOH having a concentration of 1 mol/l (corresponds to 4 % strength by weight sodium hydroxide solution, pH 14) at an operating temperature of 50 °C gives a glass surface removal of 1 mm after 25 years in a continuous flow through a glass pipeline.

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Fig. 4: Alkali attack on DURAN as a function of temperature.

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**Ecologically harmless**

Special glass is considered completely safe ecologically. Glass is manufactured from natural raw materials (sand, calcium carbonate and sodium carbonate) and does not

contain any materials which could be released and constitute a hazard to man or the environment. Glass can be reused several times and disposed of without difficulties.

**Use in microwave apparatus**

DURAN laboratory glass is suitable for use in microwave apparatus.

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**Sealing DURAN to related glasses**

With its heat expansion of  $3.3 \times 10^{-6} \text{ K}^{-1}$ , DURAN is unsuitable for direct seals with sealing metals. For seals with tungsten or tungsten glasses, the intermediate sealing glass 8448 must

be used. For seals with glasses of the molybdenum/Kovar range, two intermediate sealing glasses as shown in the figure below are recommended.


Fig. 5: Recommended technical glass seal combinations with other glasses, important metals, alloys and ceramics.

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**Optical properties**

Spectral range in which the absorption of DURAN is negligibly small: **approx. 310 to 2,200 nm**  
**absorption of brown glass up to about 500 nm**  
 SCHOTT laboratory glasses do not show any significant absorption in the visible spectral range. Thus, DURAN appears transparent and colourless. Fairly large layer thicknesses (axial view through pipes) appear greenish. When working with light-sensitive substances, the surfaces can be coloured brown.

The result is strong absorption in the short-wave region. In photochemical processes, light transmission of DURAN in the ultraviolet region is of particular importance. The degree of transmission in the UV region shows that photochemical reactions can be carried out, for example chlorinations and sulphochlorinations. The chlorine molecule absorbs in the range from 280 to 400 nm and thus serves as a carrier of the energy of radiation.


Abb. 6: Transmissionskurven für DURAN.  
 Fig. 6: Transmission curves for DURAN.

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**Physical properties and characteristics for the processing of DURAN by glassblowers**

**Maximum use temperature:** 500 °C  
**Transformation temperature DIN-ISO 3585:** 525 °C

Upon cooling from the melting temperature ( $10^1$  to  $10^3$  dPa s), the glasses go through a wide range of viscosity, which ranges from a relatively viscous region (at about  $10^4$  to  $10^8$  dPa s) to the plastic range ( $10^8$  to  $10^{12}$  dPa s) and up into the known elastic-brittle state ( $> 10^{13}$  dPa s). Thus, glass shows a constant increase in viscosity upon cooling and has no solidification point as is the case for crystalline substances. For

the processing of glass, certain sections of the viscosity range are of particular importance. In the transformation range, the elastic-brittle behaviour of the glass changes to a noticeably viscous behaviour with increasing temperature, which significantly changes all physical and chemical properties as a function of temperature. Thus, the temperature region of the transformation range determines the elimination of stress during heating and the onset of stresses during cooling of the glass. The position of the transformation range is indicated by the transformation temperature  $T_g$  DIN 52 324.

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
Fig. 7: Normal temperature dependence/viscosity curve of, for example, DURAN; viscosity ranges of important processing techniques, position of fixed points of viscosity and various limiting temperatures.

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To eliminate permanent stresses, DURAN is heated to a maximum temperature of 550 °C and kept at this temperature for 30 minutes.

For the cooling which follows, the cooling rates listed in the table below are recommended.

Temperature range	550 to 480 °C	480 to 400 °C	400 to 20 °C
Glass thickness	Cooling rate		
3 mm	12 deg/min	24 deg/min	to 480 deg/min
6 mm	3 deg/min	6 deg/min	to 120 deg/min
12 mm	0.8 deg/min	1.6 deg/min	to 32 deg/min

Table 1

These data are those of today's knowledge and are intended to provide information on our products and their possible applications. Thus, their meaning is not to ensure certain

properties of the products or their suitability for a particular application purpose. We guarantee good quality in accordance with our general conditions of sale.

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Unless stated otherwise, our products are made of DURAN borosilicate glass. Articles made of other glasses are specifically indicated:

**Important:**  
DURAN is a neutral glass of high hydrolytic resistance and belongs to glass type I, according to DAB 10, the European Pharmacopoeia and USP 23.

Name	Melt No.	Linear coefficient of expansion $\alpha$ 20/300 °C $10^{-6} \text{ K}^{-1}$	Transformation temperature °C	Density $\text{g/cm}^3$
DURAN BOROFLLOAT	8330	3.3	525	2.23
SUPREMAX	8409	4.1	745	2.57
FIOLAX	8412	4.9	565	2.34
AR-GLAS/ soda-lime-silicate glass	8350	9.1	525	2.50
SBW	8326	6.5	555	2.45

Name	Class of chemical durability		
	Water resistance DIN-ISO 719	Acid resistance DIN 12 116	Base resistance DIN 52 322/ISO 695
DURAN BOROFLLOAT	1	1	2
SUPREMAX	1	4	3
FIOLAX	1	1	2
AR-GLAS/ Soda-lime-silicate glass	3	1	2
SBW	1	1	1

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Fig. 8: Alkali attack on borosilicate glass related to temperature and calculated from the weight losses. NaOH = 1 mol/l. Exposure time: 1 h.

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Fig. 9: Alkali attack on DURAN related to pH value at 100 °C.

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## Cleaning

Laboratory glassware can be cleaned manually in a dipping bath or by machine in a laboratory dishwashing machine. A wide range of cleaners and disinfectant cleaners is available for both methods in specialist stores. As contamination during delivery of our laboratory glassware cannot be totally ruled out, may we suggest that the glass is always washed prior to use. In order to take good care of the laboratory apparatuses, they should be cleaned immediately after use at low temperature, a short residence time and low alkalinity. Laboratory apparatuses which have come into contact with infectious substances are first cleaned and then sterilised with hot air or steam. This is the only way to prevent impurities from sticking onto the glass and prevent damage of the glasses caused by any adhering chemicals.

## Manual cleaning

Wiping and scouring with a cloth or sponge, each of which are impregnated with a cleaning solution, are

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generally known. Laboratory glasses must never be treated with abrasive scouring agents or sponges, since this would damage the surface. In the dipping bath procedure, the laboratory glasses are, as a rule, placed into the cleaning solution at room temperature for 20 to 30 minutes, then rinsed with tap water and then with distilled water. Do not lengthen the time of exposure or increase the temperature except in the case of persistent impurities. In the case of laboratory glasses, lengthy times of exposure at more than 70 °C should be avoided in alkaline media, since this may result in destruction of the print.

## Machine cleaning

The machine cleaning of laboratory glasses in the laboratory dishwasher is a milder treatment than cleaning in a dipping bath. The glasses only come into contact with the cleaning solution during the relatively short washing cycles, when this solution is sprayed onto the glass surface via spray or injector nozzles.

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#### Important notes for the user!

When working with glass, the limits of this material upon change in temperature and mechanical stress must be heeded and strict safety precautions must be observed:

- Never subject glassware to abrupt changes in temperature. Thus, do not remove it from the drying oven while hot and place it on a cold or, which is worse, a wet laboratory bench. This is true in particular of thickwalled glassware, such as suction flasks or desiccators.
- Provide apparatuses with a support via a suitable stand and set up free of stress. Use, for example, PTFE corrugated bellows to compensate for stresses or vibrations.
- Glass apparatus (e. g. filtering flasks, desiccators, etc.) must be handled with care when under pressure or vacuum.
- To avoid setting up stresses in the glass, glass vessels under vacuum or pressure should not be heated from one side or with an open flame. The maximum permissible pressures quoted in the catalogue should not be exceeded.

#### Note on complaints

If in spite of our quality control measures you should have a product complaint in relation to DURAN laboratory glass, then please contact your laboratory dealer or us directly.

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- Before evacuating or pressurising glass vessels always check visually that they are not damaged in any way (heavy scratching or other mechanical damage). Damaged items should not be used under vacuum or pressure.
- Never subject glassware to abrupt changes in pressure, for example never aerate evacuated glassware abruptly. Flat-bottomed glassware (for example Erlenmeyer or flat-bottomed flasks) must not be subjected to pressure or vacuum.

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#### SCHOTT GLAS

Kundenservice Laborglas  
Hattenbergstrasse 10  
D-55122 Mainz  
Germany  
Telephone + 49-61 31-664 907  
Telefax + 49-61 31-664 05 1  
We will deal with your complaint immediately and contact you as quickly as possible.

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**Plastics used  
with laboratory glass**

Temperature resistance  
up to °C

PE	Polyethylene	- 40 to + 80
PP <sup>2</sup>	Polypropylene	- 40 <sup>4</sup> to + 140 <sup>4</sup>
PBT <sup>2</sup>	Polybutyleneteraphthalat	- 45 to + 200 <sup>3</sup>
PTFE <sup>2</sup>	Polytetrafluoroethylene	- 200 to + 260
FEP	Tetra-Fluor-Ethylene/Hexafluor-Propylene	- 200 to + 200
ETFE <sup>2</sup>	Partially crystalline ethylene/ tetrafluoroethylene copolymer	- 100 to + 180
VMQ <sup>2</sup>	Silicone rubber	- 50 to + 230
EPDM	Ethylene/propylene-diene-rubber	- 45 to + 150
PU <sup>1</sup>	Polyurethane	- 20 to + 135 <sup>4</sup>
FKM	Fluorinated rubber	- 20 to + 200
POM	Polyoxymethylene	- 40 to + 90

<sup>1</sup> All laboratory glasses coated with PU may only be cleaned while moist in order to avoid any electrostatic charge which may form.

<sup>2</sup> This plastic types are recommended by the BGA (Bundesgesundheitsamt) (Federal Department of Health) in accordance with the specifications applicable in each case.

<sup>3</sup> Changes in colour may occur at temperature stresses above 180 °C.

<sup>4</sup> For a short time without mechanical stress.

**Chemical resistance of plastics**

Classes of substances + 20 °C										
	PE	PP	PBT	PTFE/FEP	ETFE	VMQ	EPDM	PU	FKM	POM
Alcohols, aliphatic	++	++	++	++	++	+	+	++	-	+
Aldehydes	+	+	++	++	++	+		++		+
Alkaline solutions	++	++	+/-	++	++	-	++	++	-	+
Esters	+	+	+	++	++	-	++	+	-	-
Ethers	-	-	+	++	++	-	-	+	-	+
Hydrocarbons, aliphatic	-	++	++/+	++	++	-	++	++	++	+
Hydrocarbons, aromatic	-	+	++/+	++	++	-	+	++	++	+
Hydrocarbons, halogenated	-	+		++	++	-	+	-	++	+
Ketones	+	+	+/-	++	+	-	++	+	-	+
Acids, dilute or weak	++	++	++	++	++	-	++	++	++	+
Acids, trated or strong	++	++	+	++	++	-	++	+	++	-
Acids, oxidising (oxidising agents)	-	+	-	++	+	-	-	+	+	-

<p>++ = very good resistance + = good to limited resistance - = low resistance</p>
--

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## Porosity

According to their pore sizes, glass filters are divided into porosity grades from 0 to 5. Table 1 shows the porosity ranges and their main fields of application. The pore sizes indicated always refer to the largest pore in the disc. This also indicates the diameter of the particles which are only retained during filtration. Porosities are determined using the Bechhold bubblepressure method which has often been described in literature<sup>1</sup>. To achieve rapid filtration it is necessary to produce as many "passage" pores as possible without dead-ends or enclosed hollow spaces. It is in this respect that SCHOTT sintered glassware excels. The practical advantages stem from the well established properties of DURAN borosilicate glass and the special manufacturing techniques used in sintering the glass particles, which are the starting material for filter discs.



An essential condition for successfully working with glass filters in the selection of the correct porosity. Table 1 shows six porosity ranges and indicates their main fields of application.

In selecting suitable filtration apparatus it should be remembered that the nominal maximum pore size should be slightly less than the size of the smallest particles to be separated. This prevents these particles from entering the pores. It also permits highest possible flow rates without making cleaning unnecessarily difficult. This is specially important in the separation of fine-grained insoluble solid particles such as silicates and graphite.

Glass filter apparatus of porosities 3 or 4 is used almost exclusively in quantitative analysis. Different porosities are sometimes recommended for the same substances. This is explained by the fact that differing precipitation techniques for gravimetric analysis often produce different particle sizes. In case of doubt, porosity 4 is preferred as it will always allow quantitative separation of the precipitate. Porosity 3 however, has proved itself completely satisfactory in all cases for substances such as silver chloride and nickel dimethylglyoxime.

<sup>1</sup> W. Frank: GIT 11 (1967) no. 7, pages 683–688

Table 1

Porosity por.	New identification mark ISO 4793	Nominal max. pore size $\mu\text{m}$	Fields of applications, examples
	P 250	160 – 250	Gas distribution. Gas distribution in liquids at low pressure. Filtration of coarsest precipitates.
	P 160	100 – 160	Coarse filtration, filtration of coarse precipitates, gas distribution in liquids. Liquid distribution, coarse gas filtration. Extraction apparatus for coarse grain materials. Loose filter layer substrates for gelatinous precipitates.
	P 100	40 – 100	Preparative fine filtration. Preparative work with crystalline precipitates. Mercury filtration.
	P 40	16 – 40	Analytical filtration. Analytical work with medium-fine precipitates. Preparative work with fine precipitates. Filtration in cellulose chemistry, fine gas filtration. Extraction apparatus for fine-grained materials.
	P 16	10 – 16	Analytical fine filtration. Analytical work with very fine precipitates (e. g. $\text{BaSO}_4$ , $\text{Cu}_2\text{O}$ ). Preparative work with precipitates of appropriate fineness. Non-return and stop valves for mercury.
	P 1.6	1.0 – 1.6	Ultrafine filtration.

## Flow Rates

To determine possible application of glass filter discs and apparatus, it is necessary to know not only the porosity, but also the flow rates of liquids and gases. These are shown for water and air in figures 1 and 2.

The data applies to filter discs of 30 mm diameter. The flow rates for other disc diameters can be calculated by multiplying the value of a 30 mm disc by the conversion factor in table 2.

**Table 2**

Filter disc diameter mm	10	20	30	40	60	90	120	150	175
Conversion factor	0.13	0.55	1	1.5	2.5	4.3	6.8	9.7	15

### Example:

Filtration of an aqueous solution with a filter funnel of 60 mm disc diameter, porosity 4, under water-pump vacuum: Using the graph in figure 1, and with a pressure differential of about 900 mbar, read off the flow rate which is 200 ml/min. From table 2 select the conversion factor for 60 mm disc diameter and calculate the flow rate:  $200 \times 2.5 = 500$  ml/min. Since the flow rate varies greatly with the pore diameter (pore radius to the power of 4), deviations from the given value can occur. Flow obstruction can also be caused by the formation of a filter "cake" on the

sintered disc. Further flow rate variations occur when working with liquids whose viscosities differ from that of water. In these cases the flow rate is inversely proportional to the viscosity. Deviations in gas flow rates occur when filter discs are covered by a layer of water or other liquid. (Gas flow in washing processes.) More detailed information on this subject can be found in literature<sup>1</sup>.

Literature:

<sup>1</sup> W. Frank: GIT 11 (1967) no. 7, pages 683-688


Fig. 1:  
Water flow rate through filter discs of various porosities as a function of pressure differential.  
For disc diameter 30 mm.


Fig. 2:  
Air flow rate through filter discs of various porosities as a function of pressure differential.  
For disc diameter 30 mm.

## Ultrafine filtration

For ultrafine filtration filter apparatus with sintered discs of porosity 5 are used. Here the nominal value of the maximum pore size lies between 1.0 and 1.6  $\mu\text{m}$ . Experiments with bacterium prodigiosum, most commonly used for test purposes in this field, have shown that a bacterium-free filtrate is obtained with a nominal maximum pore size of 2  $\mu\text{m}$ , even when filtering very dense suspensions. A strain of nearly spherical bacteria was used. Experiments with spore-producing bacillus mesentericus led to the same result. It is interesting to note that diluted suspensions of this bacteria (15,000 to 90,000 per millilitre) could still be filtered sterile through porosity 3 filters. A bacterium-free filtrate could not be obtained however when filtering dense suspensions through them. The pores are already so narrow that all bacteria in dilute suspensions

adhere to the pore walls. In the case of dense suspensions, bacteria can still pass through once the pore walls become saturated. Real straining is only effected with a maximum pore size of 2  $\mu\text{m}$  and below; i. e. it is only here that the pores are smaller than the bacteria to be separated. Ultrafine filtration is one of the most important methods for treatment of biological solutions without using high temperatures which in many cases would lead to a change in, or decomposition of, the active ingredients in the solution. For liquid filtration, sintered glass filter funnels of standard design are used. For bacterium-free filtration of gases, e.g. in ventilation of fungal and bacterial cultures, pipeline filters are used. Here, porosity 3 is adequate, providing the space in front of the dry filter disc, on the air inlet side, is stuffed evenly and loosely with cotton wool.

## Care and Cleaning

### Temperature Change, Drying and Sterilization

Please take note of the following tips. They ensure that no inner stress occurs between the vessel and the filter discs, which can lead to the filter breaking.

1. Avoid sudden changes in temperature and unequal heating. Sintered glass filter funnels, pipeline filter tubes and other sintered glass apparatus with disc diameters exceeding 20 mm, which are to be dried or sterilized, should be placed in cold ovens or sterilizers.
2. The speed of heating and cooling should not exceed 8  $^{\circ}\text{C}/\text{min}$ .
3. Before filtering hot substances bring the filter unit slowly up to temperature in the hot cabinet.
4. Wet filter units should be heated slowly up to 80  $^{\circ}\text{C}$  and dried for one hour before increasing the temperature.

Filtration apparatus should, whenever possible, stand on its rim (stem upwards) in the oven or sterilizer. A perforated support base is advantageous since it allows air convection between the inside of the vessel and the body of the oven. If angled position of the filtration apparatus in the oven is unavoidable (pipeline filter tubes), then any point of support which is near to the filter seal position should be protected against premature heatings. This is done by using an underlay of heat-insulating material.

The apparatus should remain in the oven or sterilizer during cooling. Due to the thermal inertia of this type of oven, the cooling time is adequate.

### Cleaning new Sintered Glassware

Before using sintered glass filter apparatus for the first time, hot hydrochloric acid followed by several rinses of distilled water should be sucked through the filter disc under a good vacuum. This removes dust particles and powdered glass. It is important that each successive water

rinse be started only after the preceding one has been completely flushed through. This so-called "tear through" method must only be used for cleaning filters. It should never be adopted for preparative or analytical filtration.

## Mechanical Cleaning

**Sintered glass filters should always be cleaned immediately after use.**

If no precipitate has entered the pores, surface rinsing under the tap or with a wash bottle is often sufficient. The filter disc surface can be wiped clean with a small brush or squeegee.

(Do not remove filter medium with sharp objects.)

Where some of the precipitate has entered the pores, backflushing is necessary. In the case of porosities 0 to 2 this can be done simply by using a water tap, connecting it with rubber tubing to the stem of the sintered glass piece and allowing water

to run backwards through the filter disc. The water pressure must not exceed  $1 \text{ kp/cm}^2$ .

For porosities 3, 4 and 5 the precipitate is flushed or wiped off the disc, and water is sucked through in the opposite direction to filtration. Filters clogged by dust and dirt during gas filtration can be restored by treatment with a warm detergent solution followed by blowing through clean air from the clean side of the filter. Dirt particles are brought to the surface by the foam and removed by rising with water.

## Chemical Cleaning

If, after mechanical cleaning, some of the pores still remain clogged, or if it is desirable to make sure that no residue from previous work remains before filtering a new substance, then thorough chemical cleaning is required. The choice of solvent obviously depends on the nature of the contamination. For example:

Barium sulphate  
hot, conc. sulphuric acid ( $100^\circ\text{C}$ )  
Silver chloride  
hot ammonia liquor  
Red copper oxide  
hot hydrochloric acid  
and potassium chlorate  
Mercury residue  
hot, conc. nitric acid  
Mercury sulphide  
hot aqua regia  
Albumen  
hot ammonia liquor or  
hydrochloric acid  
Grease, oil  
carbon tetrachloride  
Other organic substances  
hot, conc. sulphuric acid with  
an addition of nitric acid,  
sodium nitrate, or potassium  
dichromate

Animal charcoal  
careful heating with a mixture of  
5 volumes of conc. sulphuric acid  
+ 1 volume conc. nitric acid to  
about  $200^\circ\text{C}$

Prolonged rinsing with water must obviously follow.

For biochemical work, cleaning with dichromate sulphuric acid should be avoided, since trivalent chromium compounds, present or newly formed by reduction, are absorbed on the surface of the filter disc. When they are released during subsequent use, biological substances can be seriously damaged. This danger is eliminated by using sulphuric acid with a nitrate or perchlorate addition. Only easily soluble reduction products are formed which can be completely removed by re-washing with water. Since hot, concentrated phosphoric acid and hot alkaline solutions attack the glass surface, they are unsuitable as cleaning agents. If they have to be filtered, and increase in pore size and thus reduced life of the apparatus is unavoidable.

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# Index

<b>A</b>		Caps	
Acid bottles, complete	32	for gas washing bottle	83
Adapters		glass	55
expansion and reduction	122	Kapsenberg	54
filter funnel	76	metal	56
multiple-, with two necks	112	plastic, screw-thread	25, 125
receiver-	121	Centrifuge tubes	45
SVS	124, 140	Ceraquick cleaning scraper	11
Adapters, stopper type	142	Claisen heads	112
Allihn, filter tubes	87	Clamps	
Allihn condensers	113	closure, porcelain	57
Analytical funnels	76	quick release	108
Apparatus		Clips for conical joint	123
desiccator, vacuum	36f.	Closure clamp, for rolled bottle	57
filtering device	83	Coil	
Aspirator bottles	33, 34	condensers, Dimroth type	114
Automatic burette	69	distillate condensers	114
		jacketed-, condensers	114
<b>B</b>		Columns	
Beakers		Vigreux	113
Berzelius	4	Condensers	
filtering (filtering jar)	4	Allihn-	113
flat flange	100	coil-, Dimroth type	114
low form	3	Coil distillate-	114
Philips	4	jacketed coil	114
tall form	3	Liebig	113
Bearings		Cones	131, 132, 135
KPG stirrer-	119	Conical flasks	5, 164
Bell jars	17	Conical joint-cones	131, 132
Bell lid, with flat flange	17	(Standardized-)	
Bends	112	standardized-joints-sockets	
Berzelius beaker	4		122, 133-135
Bloom test vessel	4	Connection pieces	120
Blow ball	69	Coupling, screwthread	
Bottles		adapter	124, 140
aspirator	33, 34	CPG stirrer bearings	118, 147
culture	54	CPG stirrer shafts	119, 148
filtering	88f.	Crucibles, filter	78
for nutrient media	54, 55	Crystallizing dishes	13
gas washing	83f.	Culture bottles	54
laboratory	19f.	Culture flasks	
reagent, DURAN	26-29	Erlenmeyer shape	55
reagent, Soda-lime	27, 29	according to Fernbach	52, 56
rolled flange	57	according to Kolle	53
roller	57	according to Roux	53
square	30, 55	penicillin	53
wash	83f.	Culture tubes	46
weighing	13	Cylinders	15, 17, 65, 120, 162-163
Breed-Demeter, square bottles	55	Cylinders, measuring	63
Buechner funnels	77	Cylinders, measuring,	
Bulb pipettes	74	blanks	162-163
Burettes		Cylinders, mixing	65
automatic	69	Cylinders, mixing, blanks	162
flasks	70		
with straight stopcock	68	<b>D</b>	
Burette stopcocks	158-159	Demeter, square bottles	55
Burette valves	157	Desiccators range	36f.
Buttress, end ground	172	Desiccators, vacuum	36
		bases	37
		plates	40
<b>C</b>		Dimroth condensers	114
Calcium chloride cylinders	120	Dishes	
Candles, filter micro	86	crystallizing	13
		evaporating	13
		Petri	51
		watch glass	14

Disposable culture tubes	47	Flat bottom	166	<b>H</b>	
Distillation thermometers	120	Flat flange	99, 172	Heads	
Distilling flasks	8	round bottom	165	Claisen-	112
Drain cock	34	evaporating	110	Drechsel-	83
Drechsel heads	83f.	evaporating, blanks	167	for washing bottles	83f.
Dropping bottles	32	filtration	75f.	still-	112
Dropping funnels	115, 116, 168	flat bottom, medium length,		Hexagonal stoppers	142, 143
Dropping funnels, blanks	168	narrow neck	5	Hose connections, plastic	127, 142
Drying tube, bent	120	flat bottom, medium length,		Hotplate	before page 175
DUROPLAN petri dishes	51	wide neck	5		
		flat bottom, ground socket	110	<b>I</b>	
<b>E</b>		Iodine determination	6	Immersion filters	86
End, Buttress, blanks	172	Kjeldahl	166	Immersion filters, micro	86
Engler distilling flasks	167	pear shape	109, 167	Intermediate receiver adapters	121
Engler flasks	8	pear shape, two necks	111	Iodine determination flasks	6
Engler flasks, blanks	167	penicillin	53		
Erlenmeyer flasks	164	round bottom, flat flange	169	<b>J</b>	
for Kapsenberg caps	54	round bottom, medium length,		Jacketed coil condensers	114
narrow neck	5	narrow neck	6, 165	Jars	15
screwthread	5	round bottom, medium length,		Jars, filtering	4
wide neck	5	wide neck	6, 165	Joint	
with ground joint	110	round bottom	6, 165	clips for conical-	123
with straight rim	54	with ground socket	110	standardized conical-,	
Evaporating flasks	110, 167	round bottom, three necks	111	cones	122
Evaporating flasks, blanks	167	round bottom, two necks	111	standardized conical-,	
Evaporating dishes	13	volumetric	63	sockets	122
Expansion adapters	141	volumetric, blanks	164		
		Flat bottom flasks		<b>K</b>	
<b>F</b>		medium length,		Kapsenberg caps	54
Feigl, spot plate	51	narrow neck	7, 166	KECK	123, 143-146
Fernbach, culture flask	52, 56	medium length,		Kjeldahl flasks	166
Filters	173-174	wide neck	7, 166	Kolle, culture flasks	53
bottles, gas washing	83	with ground socket	110	Knobbed lids	38
candles, filter, micro	86	Flat flange			
crucibles, filter	78	beakers-	100	<b>L</b>	
discs, filter	173-174	lids-	100, 170	Laboratory bottles,	
discs with glass rim	80	reaction vessels, blanks	171	pressure resistant	22
filter funnel adapter	78	Flat flanges, blanks	109	plastic coated	22
filter funnel for screwfilter	80	Flat flange lids, blanks	170-171	with DIN thread	19f.
filter candles, micro	86	Flat head stoppers	128-129	Laboratory protection plates	11
filter funnels	76	Funnels		Laboratory tube connectors	141
filter funnels, micro	87	analytical	76	Levelling bottles	32
filters, immersion	87	Buechner	77	Lids	100
filters, immersion, micro	87	dropping	115, 116	Lids, flat flange	100, 101, 170-171
filters, pipeline	81	dropping-, blanks	168	Lids, flat flange, blanks	171
filter tubes, Allihn	87	filter	75f.	Lids for filter apparatus,	
filtering bottles	88f.	filter adapter	78	Witt-type	95, 96
filtering device	75f.	filter, micro	87	Liebig condensers	113
filtering jars	4	long stem	75		
guko	91	separating	116, 117	<b>M</b>	
rubber sleeves for filter	78	separating, blanks	168	Measuring cylinders	65, 162-163
crucibles	78	short stem	75	Measuring cylinders, blanks	162-163
Filter candles	174			Measuring cylinders, low form	163
Filter discs	173-174	<b>G</b>		Measuring flasks	164
Filter disc with hole	174	Gas distribution tubes	81	Measuring pipettes	72f.
Filtration apparatus, Witt	93f.	Gas syringes	148	Metal caps	56
Flasks		Gas washing bottles	83, 162	Micro filter candles	86
culture	52f.	Glass box with cover	59	Micro filter funnels	87
distilling	8	Glass caps	55	Micro immersion filters	86
distilling, Engler	8	Ground stoppers	128, 129	Mixing cylinders	65, 162
distilling, Engler, blanks	167	Guko (rubber gaskets, conical)	91	Mixing cylinders, blanks	162
Erlenmeyer				Multi purpose cylinders	163
for Kapsenberg caps	54			Multiple adapters	
narrow neck	5			with two necks	112
screwthread	5			Museums jars	16
wide neck	5				
with ground joint	110				
with straight rim	54				

<b>N</b>		<b>S</b>		Stoppers	
Narrow neck reagent bottles	28, 29	Safety, tamper proof screw caps	31	flat-head-	128, 129, 142
Narrow/wide neck	26-29	Sample Tubes, NMR	47	flat-head-, hollow	128, 129
NMR Sample Tubes	47	Screw caps	25, 31, 135, 139	plastic	129
<b>O</b>		Screwthread adapters	140	Stoppers NS	128, 129
O-Rings	107, 172, 173	Screwthread adapter coupling	125	Straight valves	157
for desiccator	39	Screwthreads with cone	142	Supporting device for reaction	
for rubber	151	Screwthread tubes	46, 139, 142	vessels	108
flat flanges	107	for glass blowers	124	Syringes	148
Oneway stopcocks	152-155	Screwthread		<b>T</b>	
Organ storage	14	caps	25, 31, 125	Tamper-proof	31
<b>P</b>		culture tubes	46	Test tubes	43
Pear shape	109	Erlenmeyer flasks	5	Test vessel, Bloom	4
Penicillin flask	40	filtering flasks	88f.	Thermometers, Distillation-	120
Petri dishes	51, 52	laboratory bottles	19f.	Threeway stopcocks	156
Petri dishes with parts	52	tubes, with ground cone	124	Twoway stopcocks	155-156
Pipeline filters	81	Sealing	137, 140	Tube	
Pipettes, measuring	72f.	Seals	80, 125, 126	drying-, bent	120
Pipettes, bulb	74	Seals for desiccators	39	Tubes	
Phillips beakers	4	Sedimentation cones	18, 168	centrifuge	45
Plastic screw-caps	25, 125, 126, 139	Separating funnels	116-117, 168	culture	46
with aperture	125	Separating funnels, Squibb type	117	filter, Allihn	87
Plastic stoppers	129	Separating funnel blanks	168	gas distribution	81
Plastic hose connections	127, 138	Separating funnel stopcocks	159-161	NMR	47
Plastic tubulature	127	Shafts, CPG stirrer-	119	suction	87
Plate holder with standard boss	11	Side arms for stopcocks	150	test	43
Pouring rings	25	Silicone (VMQ) rubber rings		threaded with ground cone	124
Powder funnels	75		125, 140	<b>V</b>	
Pregl, micro filter funnels	87	Silicone gasket for desiccators	39	Vacuum desiccators	36f.
Pumps, water-jet	96	Slit sieves	80, 174	Valves	157
<b>Q</b>		SLK hotplates	before page 175	Valves, 90 degrees	157
Quadrupod for laboratory		Sockets	133-135	Vessels	
protection plates	11	Sockets for spirit lamps	12	Bloom test	4
Quick release clamps	108, 173	Spare pipettes		reaction, flat flange	169
<b>R</b>		for dropping bottles	32	reaction, flat flange, blanks	171, 172
Reaction vessels,		Specimen jars	16	Vigreux columns	113
flat flange	99, 169	Spherical joints	136-137	Volumetric flasks	
Reaction vessels, flat flange,		Spirit lamps	12	with ground joint	63
blanks	171	Spot plate, Feigl	51	Volumetric flasks	
Reagent bottles	26-29	Square bottle, DURAN	30	with stopper bed	164
Receiver adapters	121	Square bottle	30	<b>W</b>	
Reduction adapters	141	Square bottles, Breed-Demeter	55	Washing bottles	83f.
Releases, quick	108	Staining dish	58	Watch glass dishes	14
Rings, silicone (VMQ) rubber	107	Staining jar	58	Water connectors	137
RODAVISS	135	Staining tray	59	Water jet pumps	96
Roller bottles	57	Stainless steel handle	108	Weighing bottles	13
Rolled flange bottles	57	Standardized conical		Weighing bottles top	149
Round bottle, amber,		joints-cones	122, 131, 132	Wicks for spirit lamps	12
plastic coated	31	STERIPLAN Petri dishes	51	Wide neck reagent bottles	26, 27
Round bottom cylinders	172	Stillheads	112	Witt, filtration apparatus	93f.
Round bottom flasks		Stirrer bearings	147	Woulff bottles	108
medium length,		Stirrer bearings, CPG-	118		
narrow neck	6, 165	shafts, CPG-	119		
medium length,		Stirrer shafts	148		
wide neck	6, 165	Stopcock, single way	152-155		
medium length,		Stopcock, two way	155-156		
with ground joint	110	Stopcock, three way	156		
three necks	111	Stopcock, burette	158-159		
two necks	111	Stopcock, separating funnels			
Roux, culture flask	53		159-161		
Rubber sleeves for crucibles	78	Stopcock barrels	149		
Rubber gaskets, conical (Guko)	91	Stopcock keys	150-152		

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SCHOTT GLAS, MAINZ

# Terms of Delivery and Payment

The following terms of delivery and payment apply to all goods and services supplied by us. Should our customers' terms of purchase differ from these terms, they will only be accepted if confirmed by us in writing. Acceptance of goods supplied constitutes acceptance of these terms.

1. Prices are in German Deutschmarks at the ex works price in effect on the day of delivery. Unless agreed otherwise, invoices are to be paid immediately without deduction of any discount.

If payment is not made by the due date, legal default proceedings may be taken without further notice. We also reserve the right to charge interest at the rate charged by our own bank at the time and the total balance owing will be due for immediate payment, regardless of any other previous arrangement.

For custom-made items we reserve the right to add a reasonable price supplement and to vary the quantity supplied within reasonable limits. Any quantities in excess of the quantity ordered must be accepted. In orders for blanks and glass-to-metal seals, the amount delivered may deviate from the number actually ordered by around 10 %.

2. Every effort is made to maintain quoted delivery dates. However, due to the risks and nature of working with glass, these dates are not in any way binding, unless agreed otherwise. Our contractual obligation to our customers is subject to correct and prompt deliveries from our suppliers.

3. The point of delivery is the works where the goods are manufactured. The place for the fulfilment of payments is Mainz. The goods will be despatched through forwarding agents chosen by us. Once the forwarding agents are in receipt of the goods, the risk passes to the buyer.

4. Packing material is not returnable unless otherwise arranged.

5. Every care is taken in the execution of orders. If, in spite of this, the buyer should have any reason for complaint, notification of this must

be within 14 days of the receipt of the goods. Otherwise the goods will be deemed to have been accepted. Warranty claims for industrial equipment or parts thereof must be made within 6 months of supply in the case of those used on a single shift basis and 3 months for multi-shift use. Parts subject to wear and tear are not covered by warranty. No goods should be returned to us without our prior approval.

If, in spite of our care and attention, the goods supplied prove to be faulty, we reserve the right to either repair or replace them. In the case of damages, including consequential damages, the maximum for which we will be liable is the amount on the invoice and, in the case of glass-to-metal seals, only if more than 3 % of the batch are rejected as faulty.

6. All drawings, illustrations, dimensions and weights are approximate, unless otherwise specified. The buyer is responsible for ensuring that any working drawings he submits do not infringe the patent rights of any third party: we will not be held responsible for any resulting compensation claims.

7. Any information that we supply may not be reproduced, made available to third parties, or used for any other than the agreed purpose.

8. The buyer shall make such arrangements for storing and identification of the goods as shall be satisfactory to the seller.

Property in the goods shall not pass to the buyer until payment in full has been received by the seller for all goods sold under the contract and for all goods supplied by the seller under any other contract and until that time the buyer shall hold any goods in a fiduciary capacity on the seller's behalf.

The buyer shall not be prevented by the seller's reservation of title from selling or dealing in the ordinary course of business with any of the goods by all such sales or dealings shall be made by the buyer as agent for the seller and the buyer shall not attempt or purport to pledge or transfer the goods in any way as security to any third party until title to them has passed to the buyer in accordance with these conditions.

Any proceeds of sale of goods sold by the buyer pursuant to above conditions shall be kept by the buyer in a separate account.

The seller shall be entitled, at any time after a default has been made by the buyer on any payment by the due date, to require the buyer by written notice to place at the seller's disposal any goods to which title has been reserved by the seller under these conditions any to pay over to the seller any proceeds of sale in respect of goods which have been sold by the agent for the seller. Any action so taken by the seller shall not prejudice the rights of the seller either with respect to the goods concerned or any other goods or with respect to its action for the price.

9. In the case of sales in countries where the right of ownership described in 8 above does not have the same effect as in German law, the goods remain our property until all contractual obligations of the sale have been fulfilled, and all payments have been made. If the right of ownership does not have the same effect as in German law but the reservation of other rights over the goods is permissible, then we have the authority to exercise these rights. The buyer must be prepared to co-operate with any measures we take for the protection of our right of ownership, or the equivalent right to the goods.

10. All legal relationships formed as a result of the sales agreement are in accordance with German law. Mainz is the place of jurisdiction for all parties and all matters, including those concerning bills of exchange. If Schott appears as a plaintiff, it will also be entitled to start legal proceedings at the premises of the buyer.

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