

EG-Konformitätserklärung

EC Declaration of Conformity

Beleuchtungssysteme aus:

Illumination systems consisting of:

VisiLED MC1500 - S80-55/S80-25/S40-55/S40-10D/TLS-BF/TLS-DF/ACT/Slim Ringlight

erfüllt die Bestimmungen folgender Richtlinien des Rates der Europäischen Gemeinschaft:

- 2004/108/EG (EMV-Richtlinie) – Unter der Voraussetzung die Lichtquelle wird mit empfohlenen Transformator betrieben.
- 2006/95/EG (Niederspannungs-Richtlinie)
- 2002/95/EG (RoHS)

Die Lichtquellen sind für den Betrieb mit einer Nennspannung von 21V bis 24V ausgelegt.

Die Übereinstimmung der Lichtquellen mit den wesentlichen Schutzanforderungen obiger Richtlinien wird durch die vollständige Einhaltung folgender Normen nachgewiesen:

is in conformity with the following European Directives:

- 2004/108/EC (EMC Directive) – Provided the light source is used with the recommended transformer.
- 2006/95/EC (Low Voltage Directive)
- 2002/95/EC (RoHS)

The light sources are designed for use with a nominal voltage of 21V to 24V.

Full compliance with the standards listed below proves the conformity of the light sources with the essential protection requirements of the above mentioned EC directives:

- EN 55011:2009
- EN 61000-3-2:2006
- EN61000-3-3:1995+A1:2001+A2:2005

- EN 61326-1:2006
- EN 61010-1:2001

Nach EN 62471:2008 handelt es sich um ein Beleuchtungssystem der

Risiko-Klasse 1

Die Übereinstimmung mit den deutschen Ausführungen (DIN EN) aller oben angegebenen europäischen Normen (EN) ist ebenfalls nachgewiesen.

Die Lichtquelle trägt das CE-Konformitätskennzeichen

Acc. to EN 62471:2008 it is about an illumination system of

Risk Class 1

Full compliance with German versions (DIN EN) of all European standards (EN) listed above is also proven.

The light source bears the CE conformity mark



Mainz, November 2011

ppa.

Dr. Burkhard Danielzik

Business Unit Lighting and Imaging
Vice President Europe

i.V.

Ralf Teichmann

Business Unit Lighting and Imaging
Regional Sales Manager Europe