

# SCHOTT® Cables for Medical Industry

- One piece molded sheathing for maximum adhesion; no adhesives used.
- Autoclavable, Steris & Sterrad tested
- Medical grade A2 Fiber N.A. .64
- Ergonomic grips for easy handling with gloves.
- Private labeling available



Characteristics	
Available Lengths:	180 cm (6 feet) <b>225 cm (7.5 feet)</b> Other lengths available upon request.
Fiber Optic Bundle Diameter:	<b>4.8 mm</b> 5.0 mm Other diameter sizes available upon request.

## Universal Medical Cable Connector fits with Adapter:

ACMI

WOLF

OLYMPUS / ACMI

STORZ

OTHERS

## Specialized Customer Packaging from SCHOTT

SCHOTT can arrange for customized packaging with your company logo.

Please contact one of our sales people for more information.



SCHOTT's specialized packaging for your needs

### For more information please contact

Lighting and Imaging

**SCHOTT AG**

Otto-Schott-Strasse 2

55127 Mainz

Germany

Phone: +49 (0) 6131/66-7844

Fax: +49 (0) 6131/66-78 50

[lightingimaging@schott.com](mailto:lightingimaging@schott.com)

[www.schott.com/lightingimaging](http://www.schott.com/lightingimaging)

For your local sales contact please visit:

[www.schott.com/lightingimaging/english/contact](http://www.schott.com/lightingimaging/english/contact)

All specifications are subject to change without prior notice. This datasheet or any extracts thereof may only be used in other publications with express permission of SCHOTT.

© SCHOTT AG

**SCHOTT**  
glass made of ideas