


EasyLED Accessories

For best performance of LED illumination



Polarizing filters



Diffusor and protection window for ringlight



Color filters for spot light



Mounting adaptor

Polarization filters

- Polarization system for ringlight: turnable polarizer plus analyzer.
- Turnable polarization filter attachment for transmitted light stage.
- Turnable polarization filter with adapting thread for spot light.

Color filters for Spotlight

- Fixing at spot light with thread.
- Available colors: red, blue, green, yellow.
- Halogen filter

Accessories for Ringlight

- Diffusor with adapting thread (M110 x 1) to avoid reflections of LEDs
- Borofloat® glass protection window with adapting thread (M110 x 1).
- Protects LEDs and Microscope objective in oily and dusty environments.

Mounting adaptors for Spotlight

- Brackets for colum Ø from 25 – 35 mm
- Base for stand-alone application
- Bracket for ZEISS stand K, stand M

Accessories for EasyLED

| Polarization filters | Part. No. | Mounting |
|--|-----------|--------------------------------------|
| Polarization filter system for ringlight | 400 550 | Polarizer M110x1 Analyzer Ø 66 mm |
| Polarization filter attachment for transmitted light stage | 158 500 | Ø 84 mm |
| Polarization filter for spot light (turnable) | 600 600 | M22.5x0.5 |
| Analyzer | 158 505 | M49x0.75 |
| Analyzer | 158 510 | M52x1 |

| Accessories for Ringlight | Part. No. | Mounting |
|---------------------------|-----------|----------|
| Diffusor | 400 570 | M110 x 1 |
| Protection Window | 400 560 | M110 x 1 |

| Color filters for Spotlight | Part. No. | Mounting |
|-----------------------------|-----------|-----------|
| red | 600 650 | M22.5x0.5 |
| blue | 600 660 | |
| green | 600 670 | |
| yellow | 600 680 | |
| halogen | 600 690 | |

| Mounting Adaptors for Spotlight System | |
|--|---------|
| Bracket for column Ø 20 mm | 120 220 |
| bracket for column Ø 25 mm | 120 225 |
| bracket for column Ø 29 mm | 120 229 |
| bracket for column Ø 32 mm | 120 232 |
| bracket for column Ø 35 mm | 120 235 |
| Base for stand-alone application | 158 340 |
| Bracket for stand-alone application | 120 250 |
| Bracket for ZEISS stand K, stand M | 120 240 |