


# syriQ™ Luer Cone Syringes



Format in mm	Total Length (TL)	Glass Barrel (GL)	Inside Diameter (ID)	Outside Diameter (OD)	Flange Thickness (FT)	Flange (DD)	Flange (CD)	Flange (AD)
1.0 ml long	73.6	64.7	Ø 6.35	Ø 8.15	1.9	Ø 13.8	Ø 11.0	Ø 13.8
1.0 ml short	55.6	46.7	Ø 8.65	Ø 10.85	2.2	Ø 17.75	Ø 14.7	Ø 17.75
1.25 ml	60.4	51.5	Ø 8.65	Ø 10.85	2.2	Ø 17.75	Ø 14.7	Ø 17.75
2.25 ml	74.4	65.5	Ø 8.65	Ø 10.85	2.2	Ø 17.75	Ø 14.7	Ø 17.75
3.0 ml	92.2	83.3	Ø 8.65	Ø 10.85	2.2	Ø 17.75	Ø 14.7	Ø 17.75



tip cap



Dätwyler FM 27  
West 7025/55 grey  
West 7028/55 grey